# FORUM PanAmericano

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
International FORUM Medical Devices & Processes  
since 1999


# PROFORMANCE™ CLEANING VERIFICATION

CLEARLY VISIBLE, EASY TO INTERPRET, OBJECTIVE TESTS OF CLEANING METHODS

## LUMCHECK™

The LumCheck™ is designed as an independent check on the cleaning performance of pulse-flow lumen washers. Embedded on the stainless steel plate is a specially formulated blood soil which includes the toughest components of blood to clean.

## FLEXICHECK™

This three part kit simulates a flexible endoscope channel and is designed to challenge the cleaning efficiency of endoscope washers with channel irrigation apparatus. The kit includes a clear flexible tube, attached to a stainless steel lumen device. The test coupon is placed in the lumen and the entire device is hooked up to the irrigation port of the endoscope washer.


## HEMOCHECK™ / PROCHEK-II™

Take the guess work out of evaluating the cleanliness of instruments with the HemoCheck™ blood residue test kit and the ProChek-II protein swab test. HemoCheck™ is simple to interpret and indicates blood residue down to 0.1 $\mu$ g. The ProChek-II™ measures for residual protein on surfaces down to 0.1 $\mu$ g.

## SONOCHECK™

When the ultrasonic cleaner is supplying sufficient energy and conditions are correct, SonoCheck™ will change color. Problems such as insufficient energy, overloading, water level, improper temperature and degassing will increase the time needed for the color change. In the case of major problems the SonoCheck™ will not change color at all.

## TOSI™

Reveal the hidden areas of instruments with the TOSI™ washer test, the easy to use blood soil device that directly correlates to the cleaning challenge of surgical instruments. TOSI™ is the first device to provide a consistent, repeatable, and reliable method for evaluating the cleaning effectiveness of the automated instrument washer.

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**P**rocessing is a selection of processes to re-establish a “sterile” state for medical devices as they may come into contact with the next patient in a critical way during diagnostic or surgical procedures. How critical is processing? Apart from the specifics of the intended use, this depends on the design, materials, surface properties and the manner how we use or misuse medical devices. That is why we decide to classify – according to the Spaulding classification – some procedures as “critical”, “semi-critical” or “non-critical”. This frame, however, is debatable. If we look at endoscopy, specifically duodenoscopes, bronchoscopes and ureterorenoscopes, these may be gliding on mucous membranes but may also cut into sterile tissue.

As we discussed in the first edition of the PanAmerican FORUM, this leads to new regulations and an increased pressure on the processes executed to deliver processed medical devices. Do we clean properly while using brushes, and how? Are the diameters of the brushes correct? Do we use reusable brushes and if so, do we utilize them on one specific endoscope only or are they circulating, possibly promoting cross contamination? This has to be established in standard operating procedures (SOP) to be controlled in daily clinical routine.

The processing of medical devices is part of the intended use of reusable instruments and their accessories. Therefore, it is part of the information for use (IFU) that has to be provided by the manufacturer or distributor. Do we have these available in the processing department?

Although manual processing is most frequently used worldwide, we find an increasing number of mechanical devices (like ultrasound baths or endoscope washer-disinfectors) that support the cleaning step. We made a proposal in 2008 for consideration in the German mirror group of the standard committee (see page 31). But the discussion has only just started. This shows that processing according to the state of technology and science sometimes means more than following ISO standards because they may well be lagging behind...

Today's load carriers for WD offer defined conditions according to EN ISO 15883, *inter alia* a geometry adapted to the task and corresponding flushing. This is the basis for standardization of the performance and consumption parameters of the devices, and ultimately for the validation of automated processing.

The example of the Da Vinci system, however, shows the practical problems posed by novel medical devices for established processing departments (CSSDs). This has led to expensive consequences for many operators (upgrading processing equipment, replacement procurement, need for additional staff).

The selection and calculation of processing installations requires a detailed knowledge of the specific requirements for a particular task. This notion is given additional weight by the provisions in §5 of the updated German Medical Devices Operator Ordinance (MPBetreibV), according

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to which (processing) activities may only be carried out if the user is relying on current knowledge and the required resources are present (rooms, machinery and other work equipment). If the suitable work equipment is not available (for example a suitable WD or connectors/adapters etc.), processing must not be carried out and the medical device – while fully functional and at hand – is practically and legally not applicable.

As a result, innovative medical devices, if they are reusable, cannot be procured without specialist knowledge. The German "KRINKO 2012"<sup>\*</sup> states: "It is advisable that the medical device operator ... should also be informed about the manufacturer's requirements for the processing (according to EN ISO 17664) prior to the purchase, in order to assess the feasibility of processing and the necessary means and equipment for processing ... including the staff responsible for processing and for procurement ("QM")." Unfortunately, these guidelines are often not implemented.

Accordingly, CLEANICAL has developed a Standard Operating Procedure (SOP) for procurement of medical devices (MD) as part of our consulting work, a so-called "interface description", with the users on the one hand, and the purchasers on the other. The CSSD uses this document in advance to describe their requirements in detail, thus creating transparency regarding the department's performance spectrum.

Before the MD is procured, the MD manufacturers should be asked the following questions (according to EN ISO 17664 in the current version):

1. Has the MD successfully passed the national approval (CE mark for European countries, declaration of conformity, FDA approval in the US)?
2. Is the MD a single-use or a reusable product?
3. Is the number of processing cycles limited (documentation...)?
4. Is the manufacturer's information on processing understandable and available in the national language?
5. Is there a manual in the national language?
6. Does the manufacturer describe/offer validated procedures for processing (e.g. manual cleaning, WD cleaning, disinfection and drying processes, packaging, sterilization processes)?
7. Are there dismantling and assembly plans for MD? Are there packlists with photos? Are there going to be initial instructions for users and processors with the MD by the manufacturer?

MD are classified according to their intended use by the manufacturer and by the operator. This distinguishes medical devices from normal products. Design and processing techniques are evaluated with regard to risks that may arise (ISO 13485, 14971) and treated adequately (state of technology and science).

Processing and subsequent reprocessing must lead to safe provision of a medical device to the next patient which means correct function under hygienic conditions. Processing has manual and mechanical aspects. There is no "automated processing". Processing techniques must be appropriate. Brushes, ultrasound, and washer-disinfectors may support the processing but are not necessarily available everywhere. Therefore, the manufacturer has to provide adequate processing means for different

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settings in the information for use (IFU) following ISO 17664:2017 (soon available).

Cleaning is the removal of anorganic or organic matter. Different threshold values are discussed (e.g. 100 µg per square cm). The main problem here is the unknown recovery rate from soiled surfaces. For illustration have a look at: <http://go.beckershospitalreview.com/are-your-scopes-dirty-visual-inspecting-your-instruments>

Disinfection is the biocidal action leading to a diminution of 10<sup>5</sup> for residues. Sterilization is the biocidal action that ensures that 10<sup>6</sup> units are free of microbes. These actions are time-dependent and cannot be shortened without loss of efficiency.

Lack of time is the most “critical” for processing. If there is no time, the risks are rising for staff and patients, for the process and for the medical devices.

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# Facing the standards gap: an autoclave for the rest of us\*

J. Huijs (HEART Consultancy, Quadenoord 26871 NG Renkum, Países Bajos), E-mail: jh@heartware.nl

\* First published and presented during the congress on Appropriate health technology for low resource Settings 2014 (AHT 2014), London, organized by the Institute of Engineering and Technology, UK. [www.theiet.org/aht2014](http://www.theiet.org/aht2014)

The international standards for sterile supply are embedded in Western societies; their implementation requires a strong economy. However due to many socio-economic limitations, transferring these standards – and the resulting advanced technology – to low-income countries is bound to fail. Manufacturers are forced to comply with the standards and as a consequence appropriate sterilization equipment and related supplies for this market are hardly available anymore. A major problem is caused by the way the requirements for validation of processes have been defined. The current situation related to autoclaves (sterilizers using steam as sterilizing agent) in target health facilities was analyzed, followed by research and development and resulted in prototypes that eventually should mature to an adequate autoclave for the rest of us.

## Situation analysis of sterilization equipment in target locations

### Socio-economic context

Many health facilities have to operate under harsh conditions with poor and/or limited water supply, poor/limited or even no electricity

supply and often long distances over a very poor road network. Well trained staff is extremely scarce and companies available may be very expensive, especially for a remote hospital that may be hundreds to a thousand of kilometers away from the capital, making only the trip of the technician to come to the site already a financial challenge for the hospital.

### Mind the gap

The difference between the average spending on health care per person of a population in high income and low income countries is huge. For example in the United States this spending may be approx. USD 8200 per person per year, whereas in Malawi it is as little as approx USD 74 [2]. A difference of a factor 110! It is obvious that technology developed in a high income society assumes the infrastructure of such a society: The technology and the resulting products are embedded in their society and usually require the infrastructure and related financial resources to keep them operational. The incredibly big gap of financial resources and all its implications makes it obvious that advanced equipment developed to be used in the industrialized world is bound to fail when taken to a society where only a fraction of the financial resources are available.

### Why do mobile phones thrive and autoclaves die?

Africa is undergoing a communications revolution. Mobile phones have become accessible to major parts of the populations. Mobile phones are very high tech devices, which in spite of that are a huge success also in remote African settings, whereas advanced autoclaves are found broken down and not used in many places. A major reason here is the economy of scale and the cost per item. Mobile phones are used by the tens of thousands. For such quantities it is economically viable to set up an infrastructure to have the network and the equipment operational. Moreover the price is relatively low and thus affordable for most individuals. In addition, the telephone has a tremendous positive impact on virtually everybody's life. Autoclaves however are only used in health facilities; thus quantities are very low while the diversity of makes and models is huge. Of some make/model there may be only a single unit in the country. This makes it economically not viable for a manufacturer/supplier to set up a service network. It simply is too expensive. Thus, there is no or virtually no support for the ever more complex equipment, resulting in many broken down machines. In addition the requirements for utilities such as water and electricity has become more stringent, making the equipment even more susceptible to breakdowns.

### Standards and legislation

As for many other fields of industry and services, a wide range of standards related to sterile supply have been compiled. Main objectives are to assure safety and health of users and patients, to provide a minimum level of quality and to facilitate international trade. In Europe, the standards are issued and managed by CEN (Comité Européen de Normalisation) whereas ISO (International Standards Organization) compiles standards that are to be applied worldwide. Compliance with standards is voluntary; however meeting the requirements of standards can be used to demonstrate compliance with the relevant legislation. In the European Union this is the Medical Device Directive. In the context of the Medical Device Directive sterilizers are medical devices and thus sterilizers must comply with the directive. In order to complement the European legislation, governments may have additional national laws.


Fig. 1: Membership of ISO and CEN. No low income countries are member of the Technical Committees (TC) related to sterilization in the standards organizations and thus have no influence on the standards that are formulated

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Fig. 2, page 6: Second generation prototype of the autoclave as installed in a district hospital in Eikwe, Ghana.

Fig. 3: The camshaft controller with its 4 valves. After a stepping the handle through its distinct 12 positions, the full sterilization cycle is completed. On the right hand side the air-cooled condenser for creating the vacuum for drying

For compiling standards related to sterilization both CEN and ISO have Technical Committees (CEN: TC102, TC204; ISO: TC 198). Due to several (mostly economical) reasons, the target countries in the developing world however are not or cannot be members in these committees. Sub-Saharan Africa doesn't have a single member (Fig. 1). Therefore the specific problems related to these countries are not considered in the formulation of the standards. The content of the standards is based on the assumption that they are implemented in the member countries. The standards thus are embedded in the richer, Western societies in which the often totally different conditions of the non-member countries are not considered.

#### **Validation excludes manual control. Some observations**

It is a requirement of the standards on sterilizers that of any autoclave its processes are validated. This implies that there is to be documented evidence that the autoclave renders sterile products and that the process is reproducible. In this context human actions are not considered reproducible and thus also manual control of a machine is not reproducible. In the sequencing of a process, wherever possible,

the human factor must be excluded. However, when I, together with many other people board a bus, I put my life in the hands of the bus driver. His actions cannot be validated. We rely on him. The responsibility of the bus driver is not less than that of a person operating a sterilizer. But still busses are manually driven. And that is accepted everywhere. Why then not, at least in locations where advanced technology is not feasible, accept manual control for sterilizers by well trained operators using reliable autoclaves that result in adequate sterile goods when operated well?

State of the art of sterilization equipment: fully automatic autoclaves. Often a headache for remote hospitals

Based on the concept that human actions are not reproducible and thus that only automatic autoclaves can provide validated processes, it became a requirement that processes for autoclaves are automatic. With the current technology this implies that autoclaves are microprocessor controlled. This results in equipment with a large quantity of high-tech components with a complex control mechanism. An example: In modern standard autoclaves the operation of a valve for steam supply requires a multitude of components: the microprocessor with its software; a relay; an electric valve controlling pressurized air; and the final pneumatic valve that controls the actual supply of steam. Thus there is a chain of multiple components for a single function with each component increasing the machine's complexity and increasing the risk of breakdowns.

#### **Focus on manually operated sterilizers**

In many healthcare facilities in developing countries, the use of complex automatic sterilizers as required by the current international standards is not feasible. For many facilities they drain the already very tight hospital budget. In many cases hospitals have to resort to very basic, unsafe poorly performing manual sterilizers.

**Accepting manually controlled sterilizers**  
For this market, more robust, well-performing autoclaves are required, based on manual control and without the need of complex electronic control and related components. By allowing manual control of a sterilizer a drastic reduction of complexity and number of components can be accomplished, resulting in a more reliable essential sterile supply in the target health facilities.

#### **Approval of standards organizations and WHO/health authorities**

Allowing manual control in the context of validation would require approval of the authorities related to the publication of the standards.

It is therefore recommended to open a debate on this issue with the relevant organizations such as CEN/ISO, WHO, ICRC (Red Cross), MSF (Doctors Without Borders) and possible other stakeholders. Acceptance of manual control by such respected organizations will be crucial for the commercial success of any approved sterilizer suitable for the target market.

#### **Large potential market**

The target market for the sterilizer would at least be all low income countries. However, in large regions of middle income countries, the socio-economic situation also does not qualify for reliable operating of advanced equipment. Thus the market for such sterilizers would include a major part of the worlds' countries; which are home to a major part of the world population!

#### **Essential requirements for equipment**

For any piece of equipment, in the context of this article a sterilizer, to provide its intended function throughout its lifespan, it should meet the following requirements

##### *1. It should do its job:*

For an autoclave in our context this implies that it sterilizes materials as required in the target health facilities e. g. a common remote general district hospital or health centre in the developing world.

##### *2. The equipment meets the local socioeconomic setting where it is to be used;*

in other words, the machine should be able to run during its estimated lifespan against reasonable cost in the target location. This implies that its operation, maintenance and repairs should be compatible with the annual budget for a remote moderate hospital.

#### **Performance requirements specific to sterilizers**

For a sterilizer two fundamental performance requirements can be identified

##### *1. Processed goods must be sterile*

according to the requirements for an item to be labeled sterile as formulated in the standard EN556-1. The sterilizer thus must render products sterile inside their packaging.

##### *2. Processed goods must be dry*

as specified in standard EN 285

#### **Characteristics of goods to be sterilized in a typical district hospital**

In any health facility sterile materials are used; a part of them may be single use/disposable; a large number of items are reusable and must be cleaned, packed and resterilized after each use. In a typical hospital reusable sterile materials are used in the operating theatre, delivery room and treatment rooms. These items are to be reprocessed in the central sterile supply department (CSSD) of the hospital:

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1. Wrapped surgical instruments. They may be in drums or containers. Some surgical instruments may have lumens of limited length such as suction tubes or cannulas. Instruments for MIS (Minimally Invasive Surgery) with long narrow lumens are generally not used. Such instruments require a more advanced autoclave with a pre-vacuum process.
2. Wrapped textile goods/porous loads, such as gowns, drapes, swabs etc.
3. A combination of these

#### **Process profile**

For a steam sterilizer to render sterile products for use in a hospital environment, the process has the following phases

- 1. Air-removal.** Air acts as an insulator for heat transfer by steam. Therefore in order for the steam to be in touch with all exposed objects, all air needs to be removed from the chamber and the load
- 2. Sterilization.** The actual sterilization phase during which all micro-organisms are killed. For saturated steam the standard minimum parameters are: 134 °C for a minimum of 3 minutes; or 121 °C for a minimum of 15 minutes
- 3. Drying.** Moisture is a breeding ground for micro-organisms. That is why at the end of a sterilization cycle the load should be dry. This is accomplished by creating a vacuum, causing any residual moisture to evaporate
- 4. Air-admission.** At the end of the drying cycle, the chamber pressure is below atmospheric, thus the lid cannot be opened. Only after equalizing the pressure in the chamber to atmospheric, the lid can be opened. For preventing recontamination of the load, the air entering the chamber should pass a high quality bacteria filter.

**Ensuring adequate performance of manually operated sterilizers: research**  
Given the socio-economic context of the target countries, one of the major requirements would be that it should be possible to use manually-operated sterilizers. Automatic control can be an option, but manual control should remain possible.

It therefore was essential to get scientific evidence that autoclaves that are manually operated can actually result in sterile goods that meet the set standards for sterility and dryness. That is why already in 1999, a research programme was implemented, a collaboration of the medical engineering department of a college on technology (Hogeschool Enschede, The Netherlands) and the Dutch Institute for Public Health and the Environment RIVM [1]. The research took place in one of the laboratories of the RIVM. All procedures regarding the performance testing of the sterilizer were

done based on the prescribed performance test methods as formulated in the respective standard for such studies (EN285). During the study two types of manually operated sterilizers, at the time common in district hospitals in developing countries, were tested. Its most important conclusions:

- 1. Sterilization achieved.** The requirements for sterilization of wrapped instruments and porous loads can be met in these sterilizers, provided that adequate air-removal is achieved, by adhering to a process that was the result of the research which includes steam flushing and (three) above-atmospheric steam pulses.
- 2. Dryness achieved.** Adequate drying is feasible by using steam condensation. Therefore there is no need of an electric water ring pump, the technology that is widely used in mainstream sterilizers. Using a condenser, there is no need of electricity, there are no moving parts and only limited or no water is required.
- 3. Water ring pump eliminated.** When performing the steam penetration by above atmospheric steam pulses and creating an effective vacuum by a steam condenser, the need for an electric vacuum pump is eliminated. The use of the steam condenser reduces the vulnerability of the autoclave and reduces water consumption of the autoclave tremendously.
- 4. Manual control is complex; more simple operation is needed.** For a successful sterilization process, the operation of the 4 valves at the correct moment and correct sequence is critical. It makes operating the autoclave rather confusing and prone to operator errors. It was suggested to do further research in order to reduce the complexity for the operator.

#### **Optimizing manual process control**

The objective for this part of the research was to drastically reduce operator errors by reducing the number of controls/valves to be operated. In the conventional manually operated autoclaves on the market, each individual valve needs to be operated separately. Valves are at several locations on the machine. It thus requires operating multiple valves; usually 4 valves in the correct sequence and the correct moment. In several situations 2 valves need to be operated simultaneously. This method of operation requires a strict protocol and a lot of dedication of the staff. It thus is prone to operating mistakes. Therefore a solution was searched for mechanical control of all valves by single-knob. The research was done in collaboration with the Technical University in Eindhoven, The Netherlands in the period of 1999 until 2004. Several mechanisms were tested.

However, control by a camshaft system with multiple camdisks proved to be the most efficient, flexible, user-friendly and cost-effective solution. It is a technology that was used in autoclaves in the 1960s and 70s but was abandoned when electronic control systems took over.

*The camshaft controller has a number of advantages:*

- 1. Very simple operation:** a single knob for operating all valves, ensuring the right sequence of operation of valves which is critical for a correct process; thus chances for wrong operation are drastically reduced. Each step in the sterilization process is a distinct step of the control knob. Each disk will put each valve at its correct open or closed position. A full rotation in 12 steps results in a full, validated process for general hospital loads.
- 2. Valves can be operated mechanically by cams;** no need for electric/pneumatic valves; reducing number of components considerably. The control system can work fully independently of electricity.
- 3. Very sturdy, robust design;** built for a lifetime
- 4. Repairs** can be done by a well trained plumber
- 5. Process profiles can be changed** by changing the shape of the disks
- 6. The concept can be used for single-chamber as well as double-chamber autoclaves**
- 7. The concept of the system allows automation** by adding the required sensing and control components. An upgrade kit could be made available, including sensors, timers and a motor drive for the controller. As it would provide fully automatic control it would open the possibility to have this upgraded model pass the requirements for the current standards. And in case of a breakdown of the electronic control part of the system, it allows to resort to full manual control.

#### **Prototyping and scaling up**

5 prototypes of several generations were built and field-tested in 3 countries in Africa (Ghana, Central African Republic, Zimbabwe). They are jacketed sterilizers, equipped with the cam-disk controller and an air-cooled condenser as vacuum system. The first prototype which was installed in 2005 has been running for 8 years without repair on the control system. The next phase is to find manufacturers in order to scale up production and finally come to an adequate commercially available autoclave that the health sector in many countries is waiting for.

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## Discussion

The requirements of current standards have resulted in advanced medical equipment requiring the infrastructure and budget of an industrialized nation. There is, however, a dire need of equipment that meets the harsh local socio-economic background of many developing countries. A major cause of this situation is the lack of adequate standards and guidelines which all stakeholders involved in health care in the developing world are facing. It is crucial that this huge gap in the quest for quality, that standards organizations have committed themselves to, will be filled in. A gap that now produces sterilizers that meet all standards, but if

they should ever reach their destination in a remote hospital, they may never run a single cycle.

Through research it has been demonstrated that adequate sterilization can be performed in manually controlled sterilizers. Also by eliminating the electric water ring pump the complexity and water consumption can be reduced tremendously. A manually operated camshaft control system reduces operator errors considerably and thus improves reliability of the process. An add-on kit could make the unit fully automatic, opening the way to meet the requirements of the standards and still have the possibility to resort to manual control. The

author is in the process of forming a working group that will address the issue of an adequate sterilizer for the rest of us.

## Acknowledgements

I herewith thank all those dedicated people involved during the years of research, prototype building and field testing!

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- 1 B.Muis,ACP de Bruijn,A.W.Van Drongelen,J.F.M.M. Huijs. Optimization of the process for manually operated jacketed steam sterilizers. Central Service 2002; 10 (6): 373–384.
- 2 WHO Department of Health Statistics and Informatics (May 15, 2013) "World Health Statistics 2013".WHO, Geneva

## Rapid and Super Rapid Biological Indicators for Sterilization control

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# Water saving by use of a washer-disinfector in comparison to manual processing

Winfried Michels, Thomas W. Fengler


Fig. 2: Types of water in a washer disinfecter.

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**W**ater is the basis for all life and essential for a life with dignity. This is especially known by people who live in regions where this resource is scarce. Water is also used to carry out hygiene measures for washing, disinfection and sterilisation processes where it is important to manage the resource sparingly. Seeing this, it is astonishing that water consumption of instrument processing is not brought up more frequently.

The largest water consumption takes place

when washing surgical instruments and many assume automated cleaning to use less water than manual cleaning. Let's take a look at a small WD with capacity for 4 mesh trays: depending on the process stage we can have around 10 L inside the chamber. For the pre-wash and the main-wash, a total of two process phases, we are then using 20 L of water. To achieve the same standards of cleaning, for a manual process we would use an ultrasound bath (Michels and Roth, 2010). With an adequate capacity in a CSSD, filling the tub with detergents requires approximately 10 to 20 L of water. In comparison with a WD, ultrasound could probably use a bit less water than mechanical cleaning, but this also depends on dirtiness or a required change of the washing solution. In this process, no significant difference is observed between manual and automated processing.

This changes when rinsing off Instruments in order to remove any residual chemicals to a level which leaves no toxicological concerns for the re-use on instruments, in joints and lumens (EN ISO 15883-1). This is usually achieved safely in a WD through the two rinsing steps between cleaning and a final rinsing stage with thermal disinfection. When carrying out this rinsing process manually, devices are normally flushed with running water and joints moved shortly underwater so that the gaps and lumens are rinsed.

I rinsed 10 Crile haemostatic clamps under running water, holding and moving 2 clamps

under the stream at the same time. The used water was collected with a measuring container. For 10 instruments, a total of 3 Litres were consumed.

If we fill the above mentioned WD completely and exclusively with Crile clamps in all four mesh trays, we can fit around 160 clamps in total in a single charge. This would use around 50 Litres. For 160 clamps, the manual ultrasound will consume around 10 to 20 Litres of water. Rinsing will consume; with 10 clamps consuming 3 L, 160 consuming 3x16 L, so 48 Litres. In total, manual cleaning consumes between 60 and 70 Litres. Observing the results, we establish that processing with an automated WD requires less water.

This result favours the WD, even though these are still built based on the design of domestic washing machines from the middle of last century and having suffered no significant changes since. In places where WD were first developed, water saving was not an issue, we could draw on unlimited resources. With the right technology, the required resource consumption of these machines could be reduced significantly.

Is this everything? No. There is a hidden water consuming element in a WD, the steam condensation. If the temperature is raised above 50 °C steam will form and the volume will expand due to the temperature increase. Steam must be eliminated, but this cannot occur by simply liberating it into the work environment, as this would affect working conditions and germ-loaden aerosols could state a risk. Small WDs are equipped with a steam condenser, attached to the back of the machine, which uses a cold water spray to condense the steam and feed the waste water stream. This adds 12 to 14 Litres of water consumption to the washing process, rendering the difference between manual and automated processing very small. On very large WDs, this process is carried out with heat exchange units, creating a small vantage in water consumption. Most WDs are equipped with a drying function, where the chamber is flooded with heated air to evaporate any residual water. The resulting humid air is also treated with the steam condenser. This process has a water consumption which depends strongly on the drying time and expected final dryness of the instruments. The

Fig. 1: Left: Ultrasound bath. Right: Interior of a washer disinfector (WD), detail spray arm.

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Criteria evaluated	Manual cleaning	Automated cleaning
<i>Cleaning results</i>		
• Cleanliness of lumened instruments • Standardisation of the process • Validation of the process	difficult little questionable	easy yes possible
<i>Costs</i>		
• Procurement of machines • Cost of chemicals • Staff cost • Water consumption • Wear and tear of instruments	lower higher higher lower higher	higher lower lower higher lower
<i>Staff safety</i>		
• Risk of injury • Risk: chemicals • Risk: pathogenic microorganisms	high high high	low low low
<i>Staff time (approx.)</i>		
• Preparation of solution and soaking • Brushing • Rinsing • Drying • WD loading • WD unloading	10 25 10 10 — —	— — — — 10 5

drying process in the aforementioned WD can use up another 30 to 50 Litres of water, tilting the favour towards the manual washing process. This casts a rather disturbing image on automated WDs and the suspicion arises that they are formidable energy and water “eating” devices. Product information of these devices seems to omit information related to this specific consumption.

Surely, the installation of a WD with steam condensation is easy and fast in any areas where water, energy and waste water connections are available. However, in regards to sav-

ing water, it would be better to place the WD in a way that steam (and drying air) can be conducted to the exterior, making condensation with cold water superfluous. Only this way we can ensure that this processing method is indeed more economical than manual, not taking into account standardisation, time efficiency and reliability of automated processes. 

2. W. Michels, H. Frister, H. Pahlke, R. Frey. Testing the Cleaning Performance of Automated De-contamination Processes for Minimally Invasive Instruments (Hyg Med Volume 21/1996 – 6, pp 324–330)
3. W. Michels, K. Roth. Ringversuch zur Prüfung der Reinigung mit und ohne Ultraschall in der Praxis (Zentralsterilisation 18/2010; pp. 36–43)
4. International Standardisation Organisation (ISO). EN ISO 15883-1:2006 Washer-disinfectors – Part 1: General requirements, terms and definitions and tests (2006)

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- I. W. Michels, K. Roth, R Eibl. Assessment of cleaning efficacy based on the protein-surface relationship (Zentralsterilisation 3/2013, pp208–215)

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# The Schism in Reprocessing

Ralph Basile

Last year (2016), I attended the Association for Professionals in Infection Control and Epidemiology (APIC) conference held in Charlotte, NC [1]. APIC is always a great learning experience. The quality of the presenters and the attendees is top notch. I came away knowing more than when I arrived.

There was, however, one point of concern that has nagged me for a while and this point was again driven home during APIC – infection control practitioners, epidemiologists, and medical doctors (MDs) have what I believe is a blind spot when it comes to the role of cleaning in the overall scheme of effectively reprocessing medical devices. Certainly in presentations, speakers gave credence to the role of cleaning. But over and over again, I heard the importance of cleaning subjugated to the importance of disinfection/sterilization. No more so than in trying to meet the challenge of reprocessing flexible endoscopes.

I think this blind spot amongst these well educated, well engaged professionals comes from an institutional bias – bugs cause infections, so if you kill the bugs, there will be no infections. Many of us (and I count the readers of this publication in this group) know that cleaning is the most important stage in effectively reprocessing a device [2]. Effective cleaning has been shown to remove 99%+ of infectious agents [3], including bacteria. It “clears the way” for the final stage of reprocessing, the terminal process, to destroy what little remains in the way of microbial contaminants. Conversely, failure to effectively clean severely impairs the effectiveness of the terminal process, meaning that infectious agents, including bacteria, can survive and even thrive [4]. No clearer example of this has there been than the outbreaks associated with duodenoscopes [5].

So how does this institutional bias express itself? I heard speaker after speaker call into question the value of cleaning monitoring tests as important tools in the overall effort to render a device safe to use. What was the expressed shortcoming of these tools? They could not be directly linked to preventing patient infections [6]. Conversely, disinfection/sterilization methods were linked. But what is the criterion? When judged effective, these terminal steps destroy any remaining bacterial contaminants. Cleaning, conversely, does not. So methods for evaluating cleaning are not judged as improving patient outcomes because microbial contamination remains. But wait! If the device is not ef-

fectively cleaned, the terminal process will likely not be effective! The U.S. FDA in the latest version of the reprocessing instructions guidance document states that, “It is important to note that cleaning, disinfection, and sterilization are distinctly different processes.... Accordingly, cleaning steps should be validated separately and independently from disinfection or sterilization steps.” [7] Should not the same practice be followed by healthcare facilities?

Infection control practitioners, epidemiologists, MDs must reset their understanding of the role of cleaning. Cleaning is the first and most important step of reprocessing [8]. When assessing tools for monitoring the cleaning process, the criterion is not “did the last bug get killed”, the criterion is “was the device rendered ready for the terminal disinfection/sterilization process”? Was the device free enough of contaminants, including organic soils, to not interfere or inhibit the effectiveness of the terminal process [9]? Whether utilizing reagent tests for protein, hemoglobin, ATP, TOC, etc., the criterion is the absence of organic soils at a level that renders the device safe for handling by reprocessing personnel, reduces the chance for transfer of contaminants to the next patient and prepares it for the terminal process. A validated test method for cleaning verification testing is one that demonstrates it is capable of identifying whether a cleaned device meets the criteria [10].

During the APIC meeting, and in other venues as well, when it comes to flexible endoscope reprocessing, I often get a question like this, “What test should I use? Should I use a cleaning monitor? Or should I test for microbial contamination?” My answer is always the same, “YES!”

Do both. Each method tells you something different, and something very important. Did you get the device clean? Good first step. Did you kill the microbes? Good. Each phase of the process was successful. Why not just do the last one? Because failure to clean is sure to compromise the safety and efficacy of the device when used on the next patient. Further, buildup of residual soil over time may facilitate the formation of biofilm [11]. Microbial surveillance testing in a clinical setting is far from perfect [12]. Also, contaminants other than microbes are also a threat to the health and well being of the next patient (and healthcare workers) [13]. On the other hand, while a clean device is mostly free of microbial agents, it is not com-

pletely free, and sources other than poor cleaning can be the cause of contamination to the next patient. This includes contaminated water supply, ineffective or contaminated reprocessing equipment, or user error during the disinfection/sterilization stage [14].

Rendering a device safe and ready for use on the next patient is a complicated and complex undertaking – particularly with devices such as flexible endoscopes. To paraphrase a former U.S. Secretary of Defense, we don’t know what we don’t know [15]. Testing is the key to identifying issues in reprocessing that we don’t know we have until we discover them.

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- 1 APIC 43<sup>rd</sup> Annual Conference: <http://ac2016.site.apic.org/>
- 2 AAMI ST79: Section 7.5.1
- 3 CDC Guidelines for Disinfection and Sterilization in Healthcare Facilities, 2008: [https://www.cdc.gov/hicpac/pdf/guidelines/Disinfection\\_Nov\\_2008.pdf](https://www.cdc.gov/hicpac/pdf/guidelines/Disinfection_Nov_2008.pdf)
- 4 Alfa MJ, DeGagne P, Olson N, Puchalski T. Comparison of ion plasma, vaporized hydrogen peroxide and 100% ethylene oxide sterilizers to the 12/88 ethylene oxide gas sterilizer. *Infect. Control Hosp. Epidemiol.* 1996;17:92–100.
- 5 Epstein L, Hunter JC, Arwady MA, Tsai V, Stein L, Gribogianini M, Frias M, Guh AY, Laufer AS, Black S, Pacilli M, Moulton-Meissner H, Rasheed JK, Avillan JJ, Kitchel B, Limbagao BM, MacCannell D, Lonsway D, Noble-Wang J, Conway J, Conover C, Vernon M, Kallen AJ. 2014. New Delhi metallo-beta-lactamase-producing carbapenem-resistant *Escherichia coli* associated with exposure to duodenoscopes. *JAMA* 312:1447–1455. 10.1001/jama.2014.12720
- 6 Rutala William A. Disinfection and Sterilization: The Good, The Bad and The Ugly, Keynote Presentation at APIC 2016 43<sup>rd</sup> Annual Conference, Slide #66, <https://disinfectionandsterilization.org/files/2014/12/Slides-61-93.pptx>
- 7 FDA, Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling – Guidance for Industry and Food and Drug Administration Staff, 2015,Section V.A.
- 8 AAMI/ANSI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities, Section 7.5.1
- 9 FDA, Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling – Guidance for Industry and Food and Drug Administration Staff, 2015, Criterion 2.
- 10 Michelle Alfa. Cleaning: Recent Advances in Products and Processes and Real Time Monitoring, Disinfection, Sterilization and Antisepsis: Current Issues, New Research and New Technologies. May 2013, Edited by William A. Rutala, David J. Weber
- 11 Michelle J Alfa, Rosemarie Howie, Modeling microbial survival in buildup biofilm for complex medical devices, *BMC Infect Dis.* 2009; 9: 56.
- 12 Julia Kovaleva, Astrid Buss, Usefulness of Bacteriological Monitoring of Endoscope Reprocessing University Medical Center Groningen, Laboratory for Infectious Diseases Groningen, The Netherlands
- 13 Jahan Azizi, BS, CBET, Susan G. Anderson, MSN, MBA, CPHRM, Shawn Murphy, MS, BSN, RN, Susanne Pryce, JD, Uphill Grime: Process Improvement in Surgical Instrument Cleaning, *AORN Journal*, August 2012, Vol 96, Issue 2, 152–162.
- 14 AAMI/ANSI ST91:Flexible and semi-rigid endoscope processing in health care facilities, Introduction
- 15 Rumsfeld's Knowns and Unknowns: The History of a Quip, David A Graham, *The Atlantic*, March 27, 2014

## Lectura recomendada: Fengler/Michels “Are Processed Surgical Instruments Free of Protein?” (2001)

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# Recommended for Reading: Fengler/Michels "Are Processed Surgical Instruments Free of Protein?" (2001)

Thomas W. Fengler, Winfried Michels


Minimally invasive approaches in the 1980s brought up the need for a stepwise organization of processing, as more and more instruments proved to be problematic with regard to access of inner surfaces. Fine-mechanically sliding inserts grasped or cut tissue inside of human cavities surrounded by a hollow tube and inserted in trocars that were the keyholes for surgery. While fantastic results could be achieved for certain surgical procedures, the first being cholecystectomy, it was obvious that cleaning was difficult and residues almost certain to remain on or in the instruments. The task for processing staff became even more complex as the medical devices were capsulated and could, for a long time, not be disassembled or rinsed. We began describing the processing steps while we were managing about 1000 cholecystectomies per year at the hospital Moabit Berlin in 1994.

Some years later we wondered, how the degree of cleanliness could be measured on inner and outer surfaces of surgical instruments if the bioburden is not visible. This lead to the first clinical "multi-centre residual contamination study of processing (MRSA)", the first stage of which was performed from 1998 – 1999. The aim was to quantify amounts of proteins, which could be eluted from visibly clean surfaces of instruments used in clinical routine. The study was performed in five Central Sterile Supply Departments (CSSD) in Germany.

## The recommended article first appeared in Zentralsterilisation 9 (1) 2001, p.20–32.

Upon the 15<sup>th</sup> anniversary of the publication of the article we were invited to have it republished in an American Journal. In this context a peer review was undertaken, after which it was remarked that the organisation and the presentation of information is unorthodox, to say the least.

We believe this to be a question of perspective, really, which might result from a difference between the German and American academic style and tradition and also, possibly, from certain signs of aging of the original text. Some of the points that were emphasized 15 years ago needed to be made clear back then. Consider the importance of the process step of cleaning, as compared to the other processing steps (especially sterilization) that were in the

focus of attention. Nowadays, this may well seem a bit anachronistic. Also, the call for the execution of tests to monitor the cleaning results has been answered by the industry in the meantime (cf. Ralph Basile and Mary Ann Drosnock's article in FORUM PanAmericano No. 1). We have nevertheless decided to re-publish this article, because it summarizes the results of a study, which was not only the first study on residual contamination after automated processing (without thermal disinfection) on a large scale [1]. Furthermore, the effort has remained unrepeated for a long time, with respect to its sheer volume.

We were however not inclined to revise a 15-year-old text completely and present it as if it contained something new, because there is no new data pertaining to this study [2]. We have heeded the view expressed in the peer review criticism in so far, as we want to add a few words on the OPA method here, which was inadequately explained in the old text (in particular, there was no mention what the acronym OPA stands for; but in the bibliography): The OPA protein analysis is a quantitative method for the detection of free o-phthalaldehyde (OPA) sensitive amino groups of proteins, i.e. primary amino groups. o-phthalaldehyde or ortho-phthalaldehyde (OPA) is the chemical compound with the formula  $C_6H_4(CHO)_2$ . The molecule is a dialdehyde, consisting of two formyl (CHO) groups attached to adjacent carbon centres on a benzene ring. This pale yellow solid is a building block in the synthesis of heterocyclic com-

pounds and a reagent in the analysis of amino acids and primary amino groups of protein by forming together with a thiol component iso-indoles, which can be detected photometrically at 340 nm. The method of OPA protein analysis is described in EN ISO 15883-1 Annex C – for the validation of cleaning and disinfection processes. The OPA method reflects the current state of science and technology and is the most accurate of the methods mentioned in the EN ISO 15883-1.

The studies' results, namely that proteins are found in small amounts in the eluate of nearly every second surgical instrument, regardless of the specific instrument design, were quite sobering back in 2001. It would be gratifying, if a similar undertaking could serve to demonstrate today the progress, we, as an industry, have made in this regard.

You can access and download the article on [www.cleanical.eu/media/pdf/FreeProt.pdf](http://www.cleanical.eu/media/pdf/FreeProt.pdf)

- 1 A first, smaller attempt had been published in 1996: W. Michels, H. Fristet, H. Pahlke, R. Fery, Testing the Cleaning Performance of Automated Decontamination Processes for Minimally Invasive Instruments. *Hyg Med* 1996; 21: 324–330.
- 2 The data for our more recent examinations were obtained by analysis of validation reports, see W. Michels, K. Roth, R. Eibl, Assessment of cleaning efficiency based on the protein-surface relationship. *Zentr Steril* 2013; 21: 208–215.

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Would you like to receive the next issue of FORUM PanAmericano directly at work? Get a free copy when buying any product from your Supplier!

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Alternatively, you can send us your contact data (Medical centre, department, name and surname, direction, country and phone number) by email to: [subscriptions@cleanical.de](mailto:subscriptions@cleanical.de)

Me gustaría recibir una copia de / [I would like to receive a copy of](#) FORUM PanAmericano:

Centro médico / [Medical Center](#)

Departamento / [Department](#)

Apellidos / [Name](#)

Nombres / [First names](#)

Ocupación / [Occupation](#)

Dirección / [Address](#)

País / [Country](#)

Teléfono / [Phone number](#)

Enviar / [Send to:](#)

CLEANICAL GMBH, Augusta Hospital,  
Scharnhorststr. 3, 10115 Berlin, Alemania / [Germany](#)

Was safe. Is safe.  
Will always stay safe.

Ya era seguro. Es seguro.  
Continuará seguro.

If they are regularly inspected, they will provide faultless, enduring service for many years – without any compromises.

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# From a Catalogue of measures to a Standard Operating Procedure – Example: Processing of flexible endoscopes

Antje Hartwig, Thomas W. Fengler (CLEANICAL® GmbH, Augusta Hospital, Scharnhorststr. 3, 10115 Berlin, Germany), [www.cleanical.de](http://www.cleanical.de)

Christine Janotta (Business administrator for ambulatory medical care, ISO Auditor, QEP Multiplicator)

Thomas Richter (Care Manager at Endoscopy St. Elisabeth-Krankenhaus Leipzig GmbH, Biedermannstraße 84, 04277 Leipzig, Germany)

We present an example of a catalogue of measures for a real endoscopy clinic, which we (Hartwig/Fengler) visited as a part of our consulting activity. This catalogue is expressly not meant to be a complete depiction of the whole processing process and all significant process steps; described points are concretely those which we found that were deviated and in regard to which we expressed recommendations to the operator. Hence we are commenting methods and organisational or equipment-related deficiencies carried out concretely in this clinic. For appropriate processing, a thorough briefing by the endoscope manufacturer, as well as specialized training in processing of medical products is required. With this publication we want to demonstrate that even in a facility that follows the set of rules and operates properly, a catalogue of deficiencies can add up. For clarification: Looking at the criteria of the competent control body, ZLG handling instructions (ZLG: Central federal States Health Authority, from its German acronym), the urgency of measures that must be taken differs between

F1= “critical” (=totally unacceptable and to be rectified immediately)

F2= “severe” (=potential or immediate danger, to be rectified promptly)

F3= “others” (=not according to regulation but without immediate consequences, to be rectified in a timely manner).

However, several of these “other” objections can add up to a “severe” problem and a number of “severe” claims can create a “critical” situation.

After publication of this catalogue in Aseptica, we received, commendably, feedback in the form of further indications. We found these very helpful and do not want to withhold them from our readers. The present text constitutes a comprehensive summary of both articles “From a Catalogue of measures to a Standard Operating Procedure – Example: Processing of flexible endoscopes” (Hartwig/Fengler, aseptica Vol. 22 Issue 1, p. 16–18 [2016]) and “Endoscope-Processing: Practical Indications for implementation of the measure catalogue” (Janotta/Richter, aseptica Vol. 22 Issue 2, p. 17–21 [2016]). The further indications are presented as an Annex to our Table.

- Brushes for manual cleaning (suitable brushes: [www.mtp.de](http://www.mtp.de)): brushes with a ring on one end are not suited (exception: short valve brush)
- *Cepillos para la limpieza manual (cepillos adecuados: [www.mtp.de](http://www.mtp.de)): cepillos con un anillo en uno de los extremos no son adecuados (excepción: corto cepillo de válvula)*
- Demineralized water for processing
- *Agua desmineralizada para el procesamiento*
- Sterile water for flushing after manual disinfection
- *Agua estéril para el lavado después de la desinfección manual*
- Medical pressurized air for drying
- *Médico aire comprimido para el secado*

## Suggestions and Comments (see tables on p. 25/26)

### *Regarding no. 1: Requisites for processing: adequate rooms and facilities*

The KRINKO 2012 refers back to the RKI recommendation “Requirements for Hygiene of structural-functional design and device equipment of Endoscope facilities, recommendations by the Commission for Hospital Hygiene and Infection Prevention” (2002), in which, between others, detailed requirements for processing rooms are described. However, many CSSDs were planned and built before 2002, therefore not always compliant with regulations.

### *Pre-existing facilities*

Seldom do ambulatory and stationary facilities have access to ideal spatial conditions. Two requirements from this Regulation must always be implemented: processing must be carried out in at least one separate room and division between clean and unclean spaces must be ensured. If these prerequisites are not present, an additional room must be constructed or rented or processing must be contracted out to an external company.

Often these premises are too small, clean and unclean sides not separated sufficiently, sinks for used solutions or patient secretions not

installed, often cleaning basins are not large enough or exist in insufficient numbers. In these cases alternative solutions are necessary, i.e. strict organisational separation of processes in the clean and unclean spaces. This approach shall be described in internal operation procedures and staff is to be trained according to the special conditions. Another optimisation opportunity could be splash protection walls to reduce risk of a contamination from the unclean to clean side.

### *New facilities/planning of new facilities*

Taking into account the before mentioned recommendations is strongly encouraged. Operators (clinic owner/nursing director) are requested to include processing staff in the process.

The German health regulations require the involvement of hospital hygiene specialists. External planning companies as well as consultants require experience for the specific requirements of endoscopic processing, otherwise rooms might be planned too small or important surfaces, such as the ones for leakage tests, cleanliness or functionality controls and storage of process chemicals for daily use might be forgotten. Also, transport and other organisational solutions must be included in the planning process.

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No	Action	Procedure	Yes	No
1	Requisites for processing Adequate rooms and facilities	"Requirements for Hygiene of structural-functional design and device equipment of Endoscope facilities, recommendations by the Commission for Hospital Hygiene and Infection Prevention (KRINKO)" ( <a href="http://www.rki.de/Kommissionen/Kommission_fuer_Krankenhaushygiene_und_Infektionspravention/_Empfehlungen/Aufbereitung_von_Medizinprodukten">www.rki.de/Kommissionen/Kommission für Krankenhaushygiene und Infektionsprävention/_Empfehlungen/Aufbereitung von Medizinprodukten</a> )		
2	Requisites for Processing Adequate staff	Only medical devices (MD) that have a CE label can be used. For every MD there has to be manufacturer information for processing available. MD have to be classified before their first use according to KRINKO 2012. Flexible endoscopes are, according to KRINKO 2012, classified as "semi-critical B" and preferably processed using mechanical validated processes. If a mechanical processing is not possible, the procedure for "manual processing" from KRINKO 2012 Annex 8 should be used.		
3	Manual Processing Pre-treatment Pre-cleaning	Pre-treatment is carried out directly after examination. The endoscope is still attached to the supply unit. The insertion tube is wiped with a lint-free cloth. The distal tip is submerged into a container with cleaning solution and the suction activated in order to flush the working channel. Flushing is repeated several times. <ul style="list-style-type: none"> <li>- Alternate flushing and suction. Duration of suction/flushing: at least 20 seconds, at least 200 ml of solution.</li> <li>- After flushing/suction with cleaning solution and water, empty canal system with air until empty.</li> <li>- Turn off supply unit and separate endoscope from supply unit.</li> <li>- When depositing endoscope in transport container, deposit valves and caps with it.</li> </ul>		
4	Leakage test according to manufacturer	<b>I<sup>st</sup> Step</b> carried out before insertion into cleaning solution, because penetrating liquid damages the endoscope. The endoscope is placed securely on a dry working area and the leakage tester is attached. The pressure is increased to 160-180 mm Hg. The pressure is held during 5 minutes and can only decrease by 10mm Hg. If the pressure dips more than 10 mm Hg the leakage test has failed!  <b>In this case no processing → send to manufacturer!</b>  <b>2<sup>nd</sup> Step</b> takes place once pressure has been held constant for 5 minutes. The leakage tester remains attached to the endoscope. Now the endoscope can be inserted completely into the cleaning solution. If the pressure remains constant and no air bubbles rise up, the leakage test has been passed. Before the endoscope is separated from the tester: 1. Remove endoscope control element from cleaning solution! 2. Vent endoscope: Press button at leakage tester and release pressure. 3. Only now detach tester from endoscope. 4. Thoroughly dry connecting areas between endoscope and leakage testing device.		
5	Manual processing Washing	After a successful leaking test, insert endoscope completely and bubble-free into the cleaning solution. Fill accessible canals with cleaning solution. Add valves and distal cap. Carry out all washing steps under the surface of the cleaning solution. Using adequate brushes and lint-free disposable cloths, clean all interior and exterior surfaces thoroughly. Brush all accessible canals from proximal to distal (not back and forth, always pull the brush in the same direction). The diameter of the brush head must correspond to the diameter of the canal! Brushes with a ring on one end are not adequate! Exception: the short valve-brush. Appropriate brushes: <a href="http://www.mtp.de">www.mtp.de</a>		
6	Manual processing Rinse!!! after washing	Rinsing after manual pre-cleaning/washing: Important! Detached soil/deposits and the cleaning solution must be thoroughly removed. For this, the endoscope, as well as the cleaning brushes (reusable brushes) are put in a clean, empty basin. The detached soil must be removed from all interior and exterior surfaces, as well as any remnants of cleaning solution. Thorough rinsing of the exterior! Thorough flushing of the interior (pressurized water pistol)! Dry endoscopes well in order to avoid diluting the disinfection solution.		
7	Manual processing Disinfection	Completely submerge endoscope in disinfection solution. Fill canals completely with disinfection solution according to manufacturer specifications. The endoscope must be completely covered with disinfection solution. Remove syringe while endoscope is submerged. Cover the basin with lid. →The endoscope remains in the solution for the time indicated by chemical manufacturer specifications.		
8	Manual processing Rinsing after disinfection	After the exposure time, the endoscope shall be rinsed adequately and thoroughly with sterile-filtered water! Rinsing takes place in an empty and cleaned disinfected basin. On taps and pressured water pistols, terminal filters are installed in order to generate sterile filtered water for the final rinsing step ( <a href="http://www.aquafree.de">www.aquafree.de</a> ), see page 5. The filter life has to be taken into account. Filters must be regenerated regularly according to manufacturer instructions.		
9	Manual processing Drying	Drying is carried out from the outside with a lint-free cloth of low microbiological contamination and on the inside with medical pressurized air with max. 0.5 bar (pressurized air pistol). For example: GTI Medicare GmbH E-mail: <a href="mailto:info@gti-medicare.de">info@gti-medicare.de</a> , see page 5.		
10	Storage of processed endoscopes	Processed endoscopes are not allowed to be stored in the endoscope-suitcase! Only safely washed, disinfected and dry endoscopes can be stored. Endoscopes must be stored hanging and protected from dust, without valves. Storage time depends on the type of storage. The type of storage is to be tested microbiologically. The result of this test should then determine the storage period. Generally speaking, a processed dry hanging endoscope stored in an endoscope-cabinet can be kept there for 14 days. For a longer storage period, the instrument should be processed again before use on a patient. For storage of more than 14 days the low microbiological contamination cannot be guaranteed. The endoscope-suitcase is not a storage place. It is only used for shipment.		

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No	Action	Procedure	Yes	No
I1	Water quality feed water/ drinking water	The feed water is to be checked regularly microbiologically and chemically, at least once a year, according to drinking water regulations. Following microorganisms are not allowed to be found in the water: <ul style="list-style-type: none"> <li>- <i>Escherichia coli</i>, other enterobacteria and enterococcus</li> <li>- <i>Pseudomonas aeruginosa</i>, other pseudomonas or nonfermenters</li> <li>- <i>Staphylococcus aureus</i></li> <li>- Mykobacteria and legionella</li> <li>- Viridans streptococci</li> </ul>		
I2	Water quality sterile filtered water / final rinsing water	Final rinsing water must be tested microbiologically regularly, at least once a year. If any shortcomings are observed, we recommend quarterly microbiological controls. With good results, every six months. No microorganisms should be found in final rinsing water.		
I3	Water quality for manual processing of medical products	<b>Requirements KRINKO 2012</b> Pre-cleaning → Drinking water Rinsing → Drinking water Cleaning → Drinking water Rinsing → Drinking water Disinfection → Drinking water  After manual disinfection: Rinsing/Final rinsing only with sterile filtered water, as normally no sterilisation process follows (exception surgical use in sterile tissues ("critical"). We recommend the use of DI-water for all process steps (cleaning to disinfection). Terminal filters from Aqua free are compatible with DI-water.		
I4	Quality assurance Documentation Routine controls	For every process step of processing, a standard operating procedure (SOP) must be created: Pre-treatment/pre-cleaning, leakage testing, washing + rinsing + drying, disinfection + rinsing/final rinsing, drying, storage, preparation of cleaning solution + concentration + impregnation time → according to manufacturer information. All process steps must be documented (Who, when, what, with what, ...). Every single endoscope must be tested microbiologically at least once a year after processing → test of result quality. When deviations occur we recommend a quarterly microbiological control.		

#### **Regarding no. 2: Requisites for Processing: Adequate staff**

Staff that is entrusted with the processing of medical devices is required to show expertise in this field according to German legislation (Medical Devices Operator Ordinance, MP-BetreibV). This needs to be proven to control bodies. Another prerequisite for manual processing is the documentation, in written form, of the instructions by the endoscope manufacturer, especially in regards to design of canal systems and critical steps in processing, for example non-accessible canals.

Preferably mechanical processing should be used whenever possible. A validation of manual processing requires a great effort and does not provide the ideal current standard of endoscope processing. Disadvantages of manual processing are complex reproducibility and a potential deficient protection caused by a more intense contact. The application of manual processes (if available, mechanical processes) makes a justification for the equivalence of manual and mechanical processes necessary. A periodical check of the efficiency of the cleaning process must be carried out using appropriate test systems, e.g. ProCheck-Endo (info@procheck-shop.de).

Often business administration aspects are disregarded for manual processing, such as the effective use of staff and its cost. Manual pro-

cessing ties the staff, therefore reducing capacity of patient care and clinic organisation. What is more expensive in the end?

#### **Regarding no. 3: Manual Processing: Pre-treatment/Pre-cleaning**

Processing begins with pre-cleaning at the examination room, directly after examination. Pre-cleaning has the objective of stopping any organic, chemical or pharmacological (for example lubricant or local anaesthetics) residues from drying on the surfaces, canals or lumens of the product. Preliminary cleaning has to be carried out following manufacturer instructions, which will vary in regards to flushing quantities and the order of steps to be taken. These instructions also vary regarding the use of specific adapters for cleaning during preliminary washing, as well as which cleaning solutions to use.

#### **Regarding no. 4: Leakage Test according to manufacturer**

The leakage test takes place according to manufacturer instructions at the unclean site of the processing room. It is important that the adapter is dry when mounted on the ventilation valve. This implies a control of dryness of the adapter before every use. Many facilities dry the adapter routinely by using pressurised air, to avoid water damage on the circuit board. The endoscope must be pressurized following manufacturer instructions.

During the leakage test, the insertion tube must be stretched into different positions and bent, the bending segment also angled in every direction. Before separating the endoscope from the leakage tester and reducing the pressure, it has to be removed completely from the cleaning solution. Micro perforations are not recognized with the manual leakage test. However, an electronic leakage test can give an acoustic and visual alarm when a minimal pressure difference occurs.

It should be mentioned that a leakage tester is a hygiene trap on its own. It cannot be deposited in liquids and therefore a washing and disinfection process is not possible. Because of this, it must always be kept on the unclean side of the processing room and never in the examination room or staff lab coats.

#### **Regarding no. 5: Manual processing: washing**

Different cleaning solutions are to be applied due to their compatibility with materials and in accordance with manufacturer instructions. With endoscopes, enzymatic cleaners have proven to be effective. Annex 8, 2.2 Washing of KRINKO 2012 comments on this: "for pre-cleaning and washing of flexible endoscopes, superficially active not foaming substances (tensides), enzymatic cleaners or cleaning and washing solutions combined that

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are proven to be effective, are used." Protein fixating substances like aldehydes and alcohols must not be used in the washing process. The accounted for cleaning efficiency, to be proven by the manufacturer, is essential.

The decision whether to use a single-use or reusable brush for manual brush washing, double head brushes or brushes with a ring is to be decided by the operator. As brushes are generally more difficult to process and bents can favour damage to endoscopes, one-way brushes should be disposed of after their use. The most important criteria for a proper washing of brushes is:

- To use appropriate brushes according to manufacturer information in regards to brush-head, diameter and length;
- Do not change directions while brushing when inside a canal (do not pull and push).

Reusable brushes must be processed after every use with an ultrasound bath and subsequent disinfection. At the end of the day, after washing and disinfection, they need to be stored dry and protected from contamination. A functionality check for bents and the condition of the bristles is indispensable. Should the manual cleaning process begin right after washing of the brushes; it is recommended to remove the brush from the cleaning solution first as a loss prevention measure.

In order to fill the canal system completely with washing solution, washing adapters sourced from the manufacturer must be used. As manual processing is only possible by using said adapters it is recommended to have access to a spare adapter.

#### ***Regarding no. 6: "Manual processing: Rinse (!!!) after washing"***

To add to the previous information, it must be mentioned that for intermediate flushing (in between washing and disinfection) of all canals and cavities (including jet-canal, auxiliary flushing canal, and albaran canal) the corresponding washing adapters provided by the endoscope manufacturers should be used.

#### ***Regarding no. 7: "Manual processing: disinfection"***

For a safe disinfection of the complete canal system, the disinfection solution is required to be flushed by using the already mentioned washing adapters provided by endoscope manufacturers. If only the rinsing syringe is attached to the canal opening, this does not guarantee a safe and complete filling of the canals and therefore poses a deviation from the manufacturer instructions.

#### ***Regarding no. 8, 9, 10, 11, 12, 13: "Manual processing: Rinsing after disinfection/water quality"***

Depending on the processing step and medical device, different requirements for water quality at final rinsing are necessary. For washing and disinfection steps, water must have at least drinking quality. Additional requirements arise from different processing methods, for example DI-water in WDs. DI-water is short for "de-ionised water", also known as "de-mineralised water" or "completely desalinated water". For this, all salt is removed from the drinking water. A microbiological improvement of water quality does not occur. DI-water-cartridges or osmosis-devices are often contaminated microbiologically and therefore not adequate for a final rinsing step. For the manual final rinsing step of endoscope processing, no DI-water is required!

For the final rinsing step, in order to remove disinfection solution residues, water of at least drinking quality is used, which is free from pathogenic microorganisms. This final rinsing step has to take place in the clean part of the processing room, in a specifically prepared rinsing tub. Water used needs to be renewed after each rinsing process. According to KRINKO 2012 Annex 6, final rinsing of cystoscopes requires to be carried out using sterile or sterile-filtered water. The safest way to provide sterile final rinse water has already been named in the catalogue: sterile filters.

As described before, according to the drinking water regulations (TrinkwV), regular hygiene microbiological water tests are required. Facilities that must be controlled are building installations, water intakes of the processing facilities, drinking water storage tanks, i.e. water containers at disinfection devices, water conditioning devices, i.e. DI-cartridges, osmosis-devices and, finally, water transporting devices, i.e. spray lances or washing pistols. For a correct use of drinking water facilities, stagnation times must be excluded. This means that all water intakes must be used to obtain drinking water frequently. For pipes that are only used casually, a rinsing plan should be created or these pipes should be shut down.

#### ***Regarding no. 9 and 10: Manual processing: Drying/Storage***

To eliminate recontamination, drying is carried out on the clean side of the processing room, according to manufacturer instructions. "Quick drying of endoscopes" prior to reuse on patients should be differentiated from "drying prior to storage".

***Quick drying:*** If the MD is used on a patient shortly after processing, a "coarse" or "quick drying" process is enough.

***Drying prior to storage:*** The complete drying requires, at least a verifiable process, i.e. with special compressors of medical pressurized air that is blown through the canal system for a set time and with a set pressure (with proven drying efficiency). Compressors that are not labelled as MD with a CE-label available i.e. at construction markets do not meet the requirements.

Requirements are: pathogen free, mostly oil-free (0,1 mg/m<sup>3</sup>), is not allowed to decrease microbiological conditions of MD, pressure limited to max. 0.5 bar (manufacturer instructions).

If storage is carried out by hanging the endoscopes in dust-free storage lockers, this does not correspond to a standardised storage method. Time frame statements of 14 days should be considered as guidelines. Generally speaking, the operator is responsible for safe storage. Storage times must be determined with individual microbiological tests. A safe and possible to validate storage process can only be achieved with the use of a "dry locker" (according to DIN EN 16442:2015-5 "Storage locker with regulated environmental conditions for processes, thermos-labile endoscopes").

#### ***Regarding no 14: Quality assurance: Documentation, routine control***

To find and eliminate any weak points in processing, hygiene microbiological controls must be established in all endoscopy-units. Quality of these results must be supervised continuously. The requirement for a yearly microbiological surveillance, at least once a year for each endoscope, can be in contradiction with the goal of a validated processing step, being only a sample, which does not provide reliable evidence of processing quality.

Especially with manual processing this means that a regular control of the efficiency of processing efficiency must be carried out.

Precise test-systems as the already mentioned protein test (ProCheck-Endo) for evidence of residual contamination based on OPA-methods (SDS) is required, in accordance to the recommendations from "Guideline for validation of manual cleaning and manual chemical disinfection of MD" by DGKH, DGSV, AKI, VAH, status as of 2013. The operator is responsible for processing: he carries out a risk assessment and, based on this, control intervals, taking into account results of routine controls.

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# Load Carriers for washer disinfectors (WD)

*Working draft of the Surgical Instruments Workgroup for the ad hoc group of the DIN Working Committee AA 09 in the standards committee Medicine (standardization proposal)*

Thomas W. Fengler

## Introduction

It is often heard that standards should be regarded to be advanced legal opinions and that they represent the state of technology and science. The example of our proposal for standardization of load carriers (racks) for WD, as Part 8 of the series of standards for washer disinfectors EN ISO 15883, may illustrate the

notion that "progress is a snail", which takes a lot of time indeed.

As is well known, EN ISO 15883 is a descriptive standard designed to give manufacturers concrete advice on the description and equipment of their washer disinfectors (which is important for conformity assessment). These medical devices (as this is what they are themselves) usually contain a cleaning and disinfecting chamber, in which the medical devices to be cleaned / disinfected are exposed on differently configured load carriers in such a way, that they are flushed and rinsed through with the wash liquor as completely as possible. As can be seen from the figure, the distribution of the water in the chamber of the WD may differ dramatically, as a function of hydromechanics and the geometric arrangement of the medical devices, (see p. 13):

- Flushing water (with minimum pressure specified by the pump)
- Spray arm rinsing water (pressure-reduced via bores, depending on speed)
- Impact and reflection water (pressure reduced or low pressure)
- Drain water (low pressure)
- Atomized water (pressure-free)
- Flushing water (with branching-related reduced flushing pressure)
- Backwash water (with risk of contamination)
- Bottom water (with risk of contamination)
- No water (rinse shadow)

For validated processing, it is useful to be able to use a standardized description in order to offer the professional staff some orientation, so that clinical and sample loads are arranged correctly. But not only the terms vary (from "load carrier/rack", to "mobile unit", "(injector) trolley" or "batch cart", depending on the manufacturer) but also the geometries of the load carriers available on the market (Fig. 1).

In a letter to the chairman of the Working Committee AA 09, the author, together with the late Helmut Pahlke, had proposed to describe the different load carriers and load configurations of a WD in 2008, after various designs had been developed for surgical instruments, endoscopes, bedpans, in addition to the general requirements of the standards Part I.

This was based on our own work experience with mechanised cleaning and disinfection at the Berlin-based Surgical Instruments Workgroup. Over the past 15 years, we have developed new configurations for load carriers together with manufacturers and improved existing ones.

## Standardization and the state of the art in technology and science

On the basis of the decision of the Standards Committee NA 063-04-09 AA (N 247) of 12/2008, an ad hoc group was formed with the principal mandate to design a draft for this standardization project. Since then, not much has happened and the project has been postponed by agreement. Why?

There were two reasons for this:

1. During this time various innovative medical devices (instruments, manipulators) have been introduced to the market, which require particularly thorough rinsing with sufficient pressure, possibly in combination with a movement or removal of parts of the medical device. For those, successful cleaning has to be ensured even before for the placing on the market. The appropriate processes have to be validated for this. This takes time, as the various publications on the Da Vinci manipulation arms showed.

2. The developments of recent years, the work to be done on other parts of the standard, eg. for endoscope WD, or the debate on test soils and removal methods, made it impossible to deal with this task. There simply wasn't the time. Meanwhile, an ad hoc group is preparing a proposal that should be included in the revision of the first part, "General requirements". General requirements should of course include a clear definition of the interface between the machine and medical devices to be cleaned (with or without connectors).

It can be seen that the state of the art and science can only be partly covered by existing standards. These represent a national and/or international consensus, which in the case of the 15883 needed two ballots before gaining worldwide validity as EN ISO.

The requirements for hygiene in the processing of medical devices require a validated procedure. In EN ISO 17664: 2004 a corre-


Fig. 1: Configuration examples of load carriers for WD (upper left: Miele, right: Steelco, center: Getinge, lower left: Belimed, right: BHT)

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spondingly validated and validated procedure is already required for each automated and manual processing.

In this context, the loading of a load carrier requires knowledge about the appropriate arrangement of the medical devices to be cleaned and disinfected. The medical devices must be rinsed and flushed as completely as possible.

This is followed by thermal disinfection and final drying. Example: bowls must point upwards with the outer side, tubes must neither be too long or be bent, the flow through cavity instruments must be ensured by proper flushing and rinsing. The unambiguously defined terms of a descriptive standard help in carrying out such tasks.

In conclusion, the right way to place medical devices on a load carrier is currently not subject to standards. Nonetheless, it has to be validated. A clear description of equipment and functions helps to understand and carry out proper cleaning and disinfection.

## Cleaning Indicators for Washing and Disinfection control

Rosario - Santa Fe - Argentina | Tel. +54 341 4362015 / 4306245 / 4300616 | [www.terragene.com.ar](http://www.terragene.com.ar)

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Urs Rosenberg (Borer Chemie, Zuchwil, Suiza), E-mail: urs.rosenberg@borer.ch

Philipp Stiefel (Laboratory for Biointerfaces, Empa, Swiss Federal Laboratories for Materials Science and Technology, Suiza)

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# Evaluation of enzymatic cleaners for biofilm removal from endoscopes

Urs Rosenberg (Borer Chemie, Zuchwil, Suíça), E-mail: urs.rosenberg@borer.ch

Philipp Stiefel (Laboratory for Biointerfaces, Empa, Swiss Federal Laboratories for Materials Science and Technology, Suíça)

This work was supported by a grant of the Swiss Federal Commission for Technology and Innovation (CTI) [15429.1 PFLS-LS].

**H**ere we present some results of a project carried out for the development and evaluation of novel enzymatic cleaners targeting endoscope biofilm. In preliminary experiments we found that *Staphylococcus aureus* biofilm could be removed quite easily if the formula contained a protease. It turned out that *Pseudomonas aeruginosa* biofilm is much more resistant and hence better suited for our project. Formula optimization was carried out in a 96-well plate system by using Crystal Violet (CV) for staining of total biomass and the BacTiter-Glo assay for selective staining of viable bacterial cells. Furthermore, TOSI cleaning indicators were used to evaluate the efficacy of cleaning dried, coagulated blood. Finally, selected formulas were tested for biofilm removal from endoscope channels (PTFE tubing) as per ISO/TS 15883-5:2005 Annex F. The newly developed cleaner (deconex® PROZYME ACTIVE) containing four enzymes within a novel base formulation appeared to perform better than nine comparable commercial products.

## Introduction

Endoscopes are widely used as valuable diagnostic and therapeutic tools, however, it has been reported that medical device-associated outbreaks of infection can be linked to contaminated endoscopes more than to any other medical device. Endoscopes get in contact with different body fluids and the channels provide an ideal surface for bacterial adhesion. Viable bacterial cells have been detected in many endoscopes even after the cleaning and disinfection process. Under non-laboratory conditions, most bacteria occur in the form of biofilm. They adhere to surfaces where they are embedded in a self-produced layer of extracellular polymeric substances (EPS). These EPS provide structural integrity of biofilms and protect bacteria against adverse environmental effects such as UV irradiation, antibiotics and disinfectants, and make them much more stress-tolerant than planktonic cells. It is a huge challenge to avoid and remove biofilm, especially in moist environments like endoscope channels after use.

Endoscope associated infections often arise from incomplete removal of biofilm from sur-

faces. The long, narrow endoscope channels are difficult to clean, and the cleaning result cannot be controlled visually. Moreover, usage of harsh chemicals or high temperature could harm the sensitive materials built into endoscopes. Hence, for processing of endoscopes, mild cleaning agents are required. A possibly promising approach is to enzymatically destabilize the biofilm EPS which contains proteins, polysaccharides, lipids, extracellular DNA and other substances. Some enzymes like protease, DNase I, alginate lyase, amylase, and cellulase have been reported to support biofilm removal. Hence, including those enzymes in cleaning agents, could potentially improve the efficiency of biofilm detachment. A few enzymatic cleaners that are claimed to remove biofilm, are commercially available, but they often fail in showing expected biofilm removal efficacy in practice. One of the reasons for this failure could be the use of inappropriate test parameters during product development that might lead to overestimation of the cleaning perfor-

mance. Such test parameters include microorganism selection, biofilm formation conditions, and assay methods.

## Results & Discussion

### Comparison of different cleaners

The novel cleaner deconex® PROZYME ACTIVE, developed within this project, was compared to 9 comparable, commercially available cleaners from different manufacturers. Staining total biomass with Crystal Violet revealed that *S. aureus* was removed easily by most cleaners which contain protease, including deconex® PROZYME ACTIVE (Fig. 1) whereas the non-enzymatic cleaners (C5 & C7) were not able to remove *S. aureus* biofilm under the static conditions used. Likewise, the positive control (a mixture of 1% SDS, 1% EDTA, 1% NaOH and 0.1% NaClO), a very harsh solution, only partially removed the biofilm. In contrast, the positive control was effective against *P. aeruginosa*, but most cleaners were less efficacious in removing this biofilm (Fig. 1).


Fig. 1: Residual *S. aureus* and *P. aeruginosa* biofilm determined by Crystal Violet staining in 96-well plates after treatment with different cleaners. Commercial cleaners C1 – C9 were compared to deconex® PROZYME ACTIVE (dPA) and its corresponding base formulation without enzymes (dPA-BF). Values (optical density at 595 nm) are indicated relative to the negative control (WSH) which was set at 100%. PosC = positive control.

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A mixture of several enzymes in combination with an effective base formulation was required to get appropriate biofilm removal with a cleaner. With 90% reduction, deconex® PROZYME ACTIVE removed *P. aeruginosa* even slightly better than the positive control and similar to the best commercial cleaner C2. With 80% biofilm removal, cleaner C1 was also efficient, but the rest of the products removed less than 50% of the biofilm. To differentiate the cleaning and killing activity of the cleaners, the viability of the remaining cells after cleaning was further analyzed using the BacTiter-Glo assay. It was found that treatment with cleaners C3 and C7 exhibited a substantial biocidal effect, whereas the other cleaners did not display significant disinfecting properties (data not shown). Considering the biofilm removal ability of Cleaners C3 and C7 (Fig. 1), these cleaners rather killed bacteria instead of removing the biofilm.

The efficacy of deconex® PROZYME ACTIVE in removing artificial blood contamination (TOSI) was also compared to commercial products. Only cleaner C2 performed slightly better and cleaner C1 & C8 similar to the new formulation, while all the other commercial cleaners required longer incubation times to remove the soil. The two non-enzymatic cleaners C5 and C7 were not able to remove this artificial blood soil completely (data not shown).

#### *Performance against biofilm in endoscope channels*

Although ISO/TS 15883-5:2005 Annex F is part of a washer disinfector standard, we decided to use this method for testing deconex® PROZYME ACTIVE and other manual cleaners for their biofilm removal capabilities.

With deconex® PROZYME ACTIVE, the colony forming units (CFU) of the *P. aeruginosa* biofilm grown in PTFE tubing was reduced by more than 2 Log<sub>10</sub> (99%) compared to WSH-treated control tubing. Commercial cleaners C4 (0.28 Log<sub>10</sub> reduction) and C6 (0 Log<sub>10</sub>) did not sufficiently reduce CFU, while C1 (1.62 Log<sub>10</sub>) and C2 (1.51 Log<sub>10</sub>) were slightly less effective than the novel formulation (2.11 Log<sub>10</sub>). After treatment with cleaner C7, almost no viable bacteria were recovered (5.57 Log<sub>10</sub> reduction). The results regarding efficiency of EPS removal were similar to those for CFU, except for cleaner C7 which neither reduced proteins nor polysaccharides significantly. This indicates that cleaner C7 rather killed bacteria than removing the biofilm, which is consistent with observations in the 96-well plate experiments. Overall, the performance of the cleaners in endoscope tubing parallels that in the 96-well plate system (results summarized in Table 1).

**Table 1: Reduction relative to negative control (WSH) of biofilm biomass, viable bacteria, polysaccharides and proteins by treatment of biofilm in PTFE tubing with different cleaners.**

Reduction Determination method	Biomass OD <sub>600</sub> <sup>1</sup>	Bacteria		Polysaccha- rides Dubois <sup>4</sup>	Proteins Lowry <sup>5</sup>
		CFU <sup>2</sup>	RF <sup>3</sup>		
deconex® PROZYME ACTIVE	94.6%	99.23%	2.11	93.1%	97.9%
Cleaner C1	93.0%	97.61%	1.62	86.2%	89.4%
Cleaner C2	91.4%	96.89%	1.51	84.6%	95.1%
Cleaner C6	-10.2%	-1.32%	-0.01	-31.5%	-9.2%
Cleaner C4	19.4%	47.81%	0.28	19.2%	19.0%
Cleaner C7	-34.4%	99.99973%	5.57	13.1%	5.6%

<sup>1</sup> OD<sub>600</sub>: optical density measured at 600nm



<sup>2</sup> CFU: colony forming units

<sup>3</sup> Reduction Factor (Log<sub>10</sub> reduction compared to negative control)

<sup>4</sup> phenol-sulfuric acid method (Dubois method)

<sup>5</sup> protein quantification assay (Lowry assay)

Negative values indicate removal worse than WSH (water of standardized hardness).

**Fig. 2: a) Remaining bacteria on endoscope tubing. Nucleic acid-staining (Syto9) of *P. aeruginosa* biofilm on tubing treated with deconex® PROZYME ACTIVE (right) and with WSH (left), visualized with 20x water objective. Bacterial cells are appearing as white spots. b) Integrated density or surface coverage of residual biofilm relative to the negative control (WSH) which was set at 100%. Commercial cleaners C1, C2, C4, C6, and C7 were compared to deconex® PROZYME ACTIVE (dPA) and its corresponding base formulation without enzymes (dPA-BF).**

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Another important criterion of the standard is that cleaning should allow complete killing of all residual germs during the following disinfection step. Therefore, biofilm remaining in the tubing after cleaning, was treated with a peracetic acid-based high-level disinfectant (deconex® HLD PA / PA20). No viable bacteria were recovered from the tubing after disinfection when it was cleaned with deconex® PROZYME ACTIVE beforehand. In contrast, about 1,600 CFU per square centimeter were found in tubing treated with water (WSH) only prior to disinfection. This emphasizes the importance of an efficacious cleaning step to enable successful disinfection.

Microscopic analysis was performed to confirm biofilm cleaning. It was observed that large parts of the biofilm were removed by treatment with deconex® PROZYME ACTIVE, but not removed by WSH treatment. While the negative control displayed a dense biofilm of multiple layers (defined as 100% integrated density (ID) or surface coverage [Fig. 2b]), deconex® PROZYME ACTIVE treated sample exhibited

much lower surface coverage (17.64% ID) and almost no aggregates were observed (Fig. 2a). deconex® PROZYME ACTIVE was also found to be superior to the other cleaners regarding biofilm removal from endoscope tubing. Only cleaner C1 (44.95% ID) and C2 (31.27% ID) also displayed some bacteria-free areas, but in addition to higher coverage, more aggregates were observed (not shown). For cleaner C7 (130.35% ID), bacteria appeared blurry even though the cells were perfectly in focus. This is likely due to killing of the bacteria by destroying their membrane integrity as a similar effect was observed when disinfecting WSH-cleaned tubing with disinfectant deconex® HLD PA / PA20 (not shown).

## Conclusion

The 96-well plate assay and the ISO standard test led to matching results regarding cleaner performance. For total biofilm biomass assessment, optical density, protein and polysaccharide quantification in the ISO test correlated well with CV staining in the 96-well plate assay.

The same holds true for CFU determination compared to viable cell staining. This confirms that the 96-well plate system is an appropriate model to screen cleaner formulas for biofilm removal and for differentiation between removal and killing.

Addition of enzymes to the cleaner base formulation had a clear beneficial effect on the biofilm removal efficacy. *S. aureus* biofilm was easily removed when an active protease was present, whereas for *P. aeruginosa* single enzymes were not sufficient. An optimized enzyme mix was required to achieve efficacious *P. aeruginosa* biofilm removal. Non-enzymatic cleaners were neither effective in blood cleaning nor in biofilm removal. Among tested enzymatic endoscope cleaners, the new formula deconex® PROZYME ACTIVE showed the best biofilm removal capability plus a very good dissolution of coagulated blood. 

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# Validation of processes in the DAC Universal – combination device for cleaning, disinfection and sterilization of dental handpieces and turbines

R. Streller

## Introduction

The processing of medical devices, which are to be used in accordance with the requirements for sterile use, with suitable, validated procedures must be carried out in such a way, that the success of these procedures is comprehensively ensured and the safety and health of patients, users and third parties is not jeopardized (in Germany according to MPBetreibV).

With the validation of the processing processes, the necessary parameters are defined, which prove that the respective process (single step of the preparation, e.g. cleaning, disinfection and sterilization of medical devices) has been carried out in such a way, that the predefined specifications were definitely achieved.

## Normative background

The DAC Universal is a combination autoclave for which validation is carried out according to ISO 15883-1 and ISO 17665-1. The two standards are considered separately.

The validation of the cleaning and disinfecting process is carried out according to ISO 15883-1 (and, in Germany, according to the guideline of DGKH, DGSV and AKI for the validation and routine monitoring of automated cleaning and disinfection processes for thermostable medical devices).

The sterilization process is validated according to ISO 17665-1 or DIN SPEC 58929, taking into account EN 13060. In addition, german DGKH recommendations for validation and routine monitoring of sterilization processes with moist heat for medical devices are observed.

## Process (Fig. 1)

The process phases are defined as follows:

*Internal cleaning:* During cleaning, water is pressed through the narrow lumen of the spray canal without chemicals. This happens in pulses. *Care:* Oil is directed into the transmission channel.

*Cold external cleaning:* The cold external cleaning step with subsequent water drain is carried out three times in succession without chemicals.

*Warm external cleaning:* The warm external cleaning step with subsequent water drain is also carried out without chemicals.

*Disinfection:* During the heating phase, the temperatures inside the unit are continuously increasing in the range of 65 °C to 135 °C for at least 90 s. This corresponds to an  $A_0$  value of >> 3000 sec, which conforms to the conditions for thermal disinfection for "semi-critical B" medical devices.

*Sterilization:* After the sterilization temperature has been reached, the instruments are steam-flooded for about one minute, after which the actual sterilization takes place. This takes 3 to 18 minutes, depending on the process (134 °C/3min, 134 °C/18min or 121 °C/15min). In all processes, an  $F_0$  value of > 15 min, which corresponds to the definition of "sterile", is achieved.

After the end of the process, the lid of the DAC Universal opens, the unpacked instruments are no longer sterile.

The processing chain "cleaning, disinfection, sterilization" can be achieved successfully with every cycle of the unit.

## Validation

### Installation qualification (IQ)

The IQ is to provide and document the proof that the machine has been delivered and installed in accordance with ISO 15883-1 and its specification. Since the DAC Universal is a combination device for the cleaning, disinfection and sterilization of dental handpieces and turbines, the equipment, documentation and the supply of operating equipment must meet the requirements of ISO 17665-1 and the manufacturers' recommendations.

The installation qualification consists of:

- Check of completeness of delivery of the machine with accessories.
- The correct installation of the machine according to the installation plan.
- Setup of connections for supply and disposal
- Provision of the associated equipment.

### Operational qualification (OQ)

The OQ is to provide and document the proof that the operation of the installed equipment is within predetermined limits when used in accordance with its operating procedures. Since the DAC Universal is a combination device for cleaning, disinfecting and sterilizing dental handpieces and/or angle pieces and turbines, it is additionally proved that the installed equipment carries out the defined sterilization procedure, according to ISO 17665.

The operational qualification consists of:

- Check of the provided media such as water, air and energy.

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Internal Cleaning  
Limpieza interna

Care/Oiling  
Mantenimiento/engrasado

Source: Dentsply Sirona Brochure DAC Universal

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
**Fig. 1: The Process in the DAC Universal**

- ① Pressure test: Pressurized air is blown into the chamber for a few seconds to test for leaks.
- ② Internal cleaning: Water pulses without any process chemicals are pressed through the narrow lumens of the channels in.
- ③ Care: Oiling of the transmission mechanism
- ④ Cold external cleaning +
- ⑤ Water removal: The cold external cleaning step with subsequent water removal is performed 3 times without adding any chemicals.
- ⑥ Warm external cleaning +
- ⑦ Water removal: The warm external cleaning step with subsequent water removal is also performed without adding any chemicals.
- ⑧ Heating: The machine is heated up to sterilization temperature using steam. This part of the process is used as the thermal disinfection phase for the WD part. (The dip in the curve shows that the machine is cold and has to heat up again.)
- ⑨ Backflush: A special feature of the DAC Universal. Steam is discharged through the adapted devices. This way the instruments' lumens are flushed with steam and the air is removed. This part of the process happens before and after the actual sterilization.
- ⑩ Sterilization: The sterilizer part is normative class S – no compensation time, 180 seconds sterilization time, 3 K temperature band
- ⑪ Drying: Drying

- The correct function of the cleaning step.
- The correct function of the disinfection step.
- The correct function of the sterilization step.

#### Performance qualification (PQ)

The PQ is to provide and document the proof that the device, as installed and operating according to the operating procedures, operates consistently in accordance with the predetermined criteria, thereby obtaining products which meet their specifications. In addition, it must be demonstrated that the product has been exposed to the specified sterilization process.

In order to demonstrate the reproducibility, one program has to undergo three test runs and all other programs must be tested once.

The following parameters are relevant for performance qualification:

1. Cleaning performance, ISO 15883-1/-2, ISO / TS 15883-5
2. Disinfecting performance, ISO 15883-1/-2
3. Sterilization performance, ISO 17665-1, DIN SPEC 58929

#### Testing the cleaning performance

The cleaning performance test is carried out with a defined test soil according to ISO/TS 15883-5 and real contaminated medical devices. Because of the small chamber volume, the number of test pieces is reduced to three per process.

To evaluate the cleaning performance, the process has to be interrupted after the last cleaning step and the test instruments have to be taken out, according to ISO 15883-1/-2.

#### Testing the disinfection and sterilization performance

For the thermal evaluation, ISO 15883-1/-2 and ISO 17665-1 must be considered.

Thermal disinfection is defined in ISO 15883-1/-2. The classification of the contra-angle pieces as semicritical B requires an A0 value of 3000sec. Sterilization is assessed according to ISO 17665-1. Each program must be considered individually (134 °C/3min, 134 °C/18min, 121 °C/15min). An F<sub>0</sub> value of 15 min is reliably achieved in each of these processes.

Disinfection and sterilization take place, one after the other, in the same process. However, the two temperature profiles, disinfection (Fig. 1 ⑦) and sterilization (Fig. 1 ⑨) can be recorded in one cycle and then evaluated separately. The process is fully recorded for thermal evaluation.

#### Validation equipment

A simple insertion of the temperature sensors or the test soils is not possible in the DAC Universal because of the special adaptation to the handpieces and contra-angle handpieces. In addition, instruments and any auxiliary means must be very well attached due to the mechanical load during cleaning.

On the lid of the DAC Universal are the adapters for the handpieces and contra-angle handpieces and the supply lines of these adapters. Thanks to the new AL 3300 system from ebro it is now possible to carry out the thermal test (disinfection, sterilization) and the cleaning performance test during validation with the AL 3301 on the standard cover used by the user.

The aim of this system is performance qualification with the equipment parts used in real processes without modifications to the DAC Universal.

#### Performance qualification of the cleaning performance according to ISO 15883-1, ISO/TS 15883-5

Adapter AL 3301 from ebro with test soil by SMP (Fig. 3). The test soil and the real contaminated handpieces and contra-angle pieces can be tested simultaneously in this system (Fig. 2). Due to the special adaptation of the AL 3301 to the laboratory-supplied reproducible test contamination, the

**Fig. 2**

**Fig. 4**

**Fig. 3**

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© Redactor jefe y mhp-Verlag GmbH 2017

All rights reserved. No part of this publication may be reproduced, transmitted in any form or by any means without the prior written permission of mhp-Verlag GmbH. One photocopy of individual articles or parts thereof may be made for personal use only. © Editor and mhp-Verlag GmbH 2016

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The use of general descriptive names, trade names, trademarks, etc., even if these are not especially identified, is not to be taken as sign that such names may be used freely by anyone. Editor in Chief and publisher do not guarantee for articles published outside the editorial part of the journal (esp. advertising, industrial information). Authors who are mentioned with their full name and who are not member of the editorial staff publish articles in their own responsibility.

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external and internal cleaning can be connected to the same adapter. This leaves slots for the real contaminated handpieces and contra-angle pieces.

The complicated method of soiling with blood, which is inconvenient in a dental practice setting, is no longer necessary. The validation of WD has long been based on contaminated Crile clamps. The test soil offered by SMP (Fig. 3) consists of a cannula similar to a spray channel and a test specimen, which corresponds to the surface of a handpiece and contra-angle piece.

#### ***Performance qualification of disinfecting performance according to ISO 15883-1/-2 and sterilization performance according to ISO 17665-1***

For the thermal test, the AL 3300 set is required for the data loggers from the EBI 11 series SL 3300 validation set (Fig. 4).

This ensures reliable, reproducible positioning of the EBI 11 temperature sensors.

The thermal load in the form of M12 x 100 stainless steel screws, conforming to EN 13060, is included in the AL 3300 and replaces the instruments in this test. The mass of the screws is higher than that of the handpieces and contra-angle pieces and can thus be regarded as a „worst case“ load (Fig. 5).

The thermal disinfecting performance and the sterilization performance take place in the same process one after the other and can be recorded in one cycle. However, the two process parts have to be evaluated separately (Fig. 1).

With the PCD testpiece, which is contained in the AL 3300, and also has a temperature sensor in position, it is checked whether the sterilization performance within the instruments is sufficient.

#### **Conclusion**

The new AL 3300 adapter set from ebro for the validation of thermal processes in the DAC Universal can be used to validate the processes. Reproducibility is provided by aid of defined positions of the data loggers and their sensors and is further increased by the defined test soil together with the AL 3301. The test with the standard cover results in a validation under real conditions.

By using the adapter set AL 3301, the test soil can be tested together with real contaminated instruments in one cycle. This means that test cycles and, through the use of test specimens, preparation time can be saved. 

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# CSSD Process Technician – an approach to international education

Johannes Guidé

The KLS Martin Group exports its products to more than 100 countries worldwide, with the focus outside the EU. As an expert in the areas of surgical instruments, sterilization containers and instruments processing, the CSSD/RUMED is also a central point of contact for the employees of the KLS Martin Group. Above all in the international environment, the processes and the education levels of the skilled workers are often not comparable with those in Germany or Europe. For this reason, the KLS Martin Group decided to set up a new educational strategy.

## CSSD Process Technician

A key stone of this strategy is the 14-day CSSD Process Technician Course at the KLS Martin World in Tuttlingen, Germany. The course itself is certified by the DGSV e.V. (German Society for hospital sterile supply) as part of the foreign language courses. In addition, the American Society IAHCSSM also contributes continual education points, which can be registered there, if the course is successfully completed.

The course itself consists of theoretical and practical approaches to convey the knowledge

in various forms to the participants. The basis for this is the curriculum of the DGSV e.V. with a final practical and theoretical examination. The combination of speakers from active practice and industry enables extensive expert knowledge being transmitted to the participants.

The concept was first introduced in English in September 2016. The 20 participants from a total of 12 countries successfully completed the course and evaluated the course as an excellent concept for further training and professional exchange.

## Future prospects

The success of the first course shows that this training concept meets the wishes and needs of the users. For this reason, the KLS Martin Group offers the same course content from 09.-20.10.2017 directly after the World Congress of the WFHSS/DGSV e.V. in Bonn. In addition, this concept is offered as a one-week basic course starting in 2017 in various countries, including Chile, China and Australia. The aim is to enable an even better local training in national language.

For further information [www.klsmartin.com/events](http://www.klsmartin.com/events), by scanning the side placed QR code or you can simply register by using our new app for Android or IOS, which can be downloaded using the QR Code below.

KLS Martin Surgical App  
Surgical Instruments, Education, Reprocessing

**kls martin**  
GROUP

# Last but not Least: Interim Report “ReSt” – residual contamination in endoscope channels

Thomas W. Fengler, CLEANICAL Investigation & Application

## Introduction:

To date there are hardly any reliable data as to the existence and extent of the residual contamination in the working channels of flexible endoscopes. In cooperation with two endoscope manufacturers and SYNLAB testing laboratory, CLEANICAL Berlin has been involved since April 2016 in a long-term study on the residual contamination of flexible endoscopes. Here, the endoscope manufacturers thankfully provide defective endoscopes (after real use) for a destructive investigation. The focus is on the health protection of the employees working in the repair departments of the manufacturers.

We would like to thank Mrs. A. Schmidt-Warnecke, SYNLAB, Berlin.

## Subject:

What degree of residual contamination can be found inside the working channels of clinically used flexible endoscopes? What does this mean for the manufacturer/repair company?

## Material:

Gastroscopes and ureteroscopes from two manufacturers

## Method:

Destructive inspection of all work channels after at least 4 weeks of storage; microbiological

examination: comminution and elution of the work channels in NaCl; filtration of the eluate and application to various selective media.

## Preliminary results:

Frequent findings of apathogenic ubiquitous microorganisms; few finds of potential pathogens (see table).

A further development of the study design, as well as an analysis of the results and possible recommendations for the protection of health will take place in a later phase, when more sample results are present (currently n = 27). Microbiological investigations have to take into account the circumstances and reduce costs as far as possible. In doing so, comparisons and a

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weighing of goods are carried out (professional cleaning, disinfection necessity). In the same way, an incidence of microorganisms on the instruments must be considered in relation to the other (the normal, everyday) microbial load. Our own study outlined here ("ReSt" – «Returns stripping» or «Residual Contamination Study») is merely an effort in diligence. Results should be evaluated with the necessary caution. We owe it to the safety of the employees to make that effort! However, this does not automatically mean, that additional measures (e.g. personal protective equipment) are necessary.

Regardless of this, the nature and efficiency of the processing of flexible endoscopes must be discussed (see the article of Basile and Drosnick in FORUM PanAmericano No.1). The first is done with regard to employee safety, the latter is for patient safety! And there, at the site of use, destructive examinations of the endoscopes are not possible; it is necessary to improve the recovery rate of residues from the working channels.

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**Table 1: Results of the microbiological examination. The sampling locations are in each case 4 or 5 regions on a total of 6 defective endoscopes after actual use on the patient (e.g. working channel, severed distal working end, light guide)**

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# “Cleanliness at Work” –

with quality products from **KARL STORZ**

More information on  
“Hygienic OR Products”  
from KARL STORZ at  
[www.karlstorz.com](http://www.karlstorz.com)

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
## Quality cycle of instrument reprocessing

→ ISO 17664

**cleanical® gmbh**

**LCMP** GmbH  
Medical Material Processing

**STÖRZ** **Miele**  
KARL STÖRZ – ENDOSKOPE  
PROFESSIONAL

**borer**  
reinforced cleaning solutions

**TERRAGENE®**

**W healthmark**

**kls martin**  
GROUP

**-ebro-**  
a xylem brand